


Method Development Guidelines

Primary structure analysis methods

This section on column selection strategy for primary structure analysis provides some critical details on method development for mAb, proteins, and peptides.



SPP = superficially porous particle, TPP = totally porous particle

Start at low pH with simple aqueous/organic gradient

Typically, a water:acetonitrile with 0.1% trifluoroacetic acid (TFA) gradient is used to elute all components of interest. A typical high resolution gradient on a 300 Å pore size column requires 30 to 50 min. An AdvanceBio RP-mAb column requires a shorter analysis time and a higher flow rate, and still provides exceptional resolution. To improve resolution, increase the gradient time, decrease column length, or increase flow rate. For LC/MS methods, TFA can reduce detector sensitivity and is often replaced with ammonium formate/formic acid.

Optimize sample solubility

For best peak shape and recovery at any pH, it is important to completely solubilize a sample. Highly acidic or neutral solvents can be used with AdvanceBio RP-mAb, ZORBAX 300 Å StableBond, Poroshell 300 StableBond, and AdvanceBio Peptide Mapping, while neutral solvents and dilute bases can be used with ZORBAX 300Extend-C18 and Poroshell 300Extend-C18.

Solvent choices to solubilize proteins and peptides

Water/phosphate buffer

Dilute acid (TFA, acetic acid or HCl)

Neutral pH, 6–8 M guanidine-HCl or isothiocyanate

Acetic acid 5%/6 M urea

Dilute acid + aqueous/organic solvents (ACE, MeOH, THF)

Dilute base (ammonium hydroxide)

DMSO or 0.1%–1% in DMSO

Formamide

Weakest

Strongest

Increase the temperature

Separations of proteins and peptides are influenced by temperature and higher column temperature can dramatically improve both resolution and recovery of proteins and hydrophobic and aggregating peptides.

AdvanceBio RP-mAb: Up to 90 °C
ZORBAX 300 StableBond, Poroshell 300 StableBond: Up to 80 °C
AdvanceBio Peptide Mapping: Up to 60 °C

Optimize mobile phase pH

Try mid and high pH if low pH does not work

If an optimized, low pH method does not provide an ideal separation, then mid or high pH mobile phase can be used. At high pH, selectivity is often very different because acidic amino acids become negatively charged and some basic amino acids may lose their charge. ZORBAX 300Extend-C18 is an excellent choice for mid to high pH separation.

Column:	ZORBAX 300Extend-C18	Gradient:	5–60% B in 30 min
	4.6 x 150 mm, 5 µm	Temperature:	25–30 °C (<60 °C)
	773995-902	Flow rate:	1 mL/min
Mobile phase:	A: 20 mM nH ₄ OH in H ₂ O		
	B: 20 mM nH ₄ OH in 80% ACN		