

# Structural Analysis of Labeled *N*-Glycans from Proteins by LC-MS/MS Separated Using a Novel Mixed-Mode Stationary Phase

Udayanath Aich<sup>1</sup>, Julian Saba<sup>2</sup>, Xiaodong Liu<sup>1</sup>, Srinivasa Rao<sup>1</sup>, and Chris Pohl<sup>1</sup>

<sup>1</sup>Thermo Fisher Scientific, Sunnyvale, CA, USA; <sup>2</sup>Thermo Fisher Scientific, Mississauga, ON, Canada

## Key Words

GlycanPac AXH-1, LC/MS, LC-MS/MS, HILIC, WAX, mixed-mode, labeled *N*-glycans, UHPLC, MS detection, Q Exactive, charge, SimGlycan software

## Abstract

This application note describes the liquid chromatography-mass spectrometry (LC/MS) analysis of fluorescently labeled *N*-glycans released from proteins. The chromatographic separation is carried out with a novel Thermo Scientific™ GlycanPac™ AXH-1 (1.9 μm, 2.1 × 150 mm) column for high-resolution and high-throughput analysis of glycans. This column possesses unique selectivity that provides separation of glycans based on charge, size, and polarity. MS and MS/MS analyses are performed using a Thermo Scientific™ Q Exactive™ hybrid quadrupole-Orbitrap™ mass spectrometer in negative ion mode to provide detailed structural information of *N*-glycans released from proteins.

## Introduction

Glycans are involved in a wide range of biological and physiological processes, including recognition and regulatory functions, cellular communication, gene expression, cellular immunity, growth, and development [1]. Glycans are commonly investigated as important species in therapeutic protein drug development because there is strong evidence that bioactivity and efficacy are affected by glycosylation [2]. Commonly, both the structure and types of glycans attached to the proteins are examined. Understanding, measuring, and controlling glycosylation in glycoprotein-based drugs, glycan content of glycoprotein products, as well as thorough characterization of biosimilars has become increasingly important.

The structures of glycans are highly diverse, complex, and heterogeneous due to post-translational modifications. This makes it challenging to comprehensively characterize glycan profiles and determine their structures [3]. It is therefore essential to separate all isomeric, charge, and branching glycan variations to understand the detailed structure of the glycans by LC-MS/MS methods.

Various HPLC separation modes have been used for the analysis of glycans, including normal phase (NP) or hydrophilic interaction (HILIC) chromatography, ion-exchange (IEX) chromatography, and reversed-phase (RP) chromatography. Because they are highly



hydrophilic, polar substances, neutral glycans are commonly separated using amide HILIC columns, such as the Thermo Scientific™ Accucore™ 150-Amide-HILIC column [4], which separates glycans by hydrogen bonding, resulting in a size and composition-based separation. This type of column is particularly useful for the separation of glycans released from monoclonal antibodies, of which the majority are neutral [5].

Based on novel mixed-mode surface chemistry, the GlycanPac AXH-1 column combines both weak anion-exchange (WAX) and HILIC retention mechanisms for optimal selectivity and high resolving power [6]. The WAX functionality provides retention and selectivity for negatively charged glycans, while the HILIC mode facilitates the separation of glycans of the same charge according to their polarity and size. As a result, the GlycanPac AXH-1 column

**Thermo**  
SCIENTIFIC

provides unparalleled capabilities for glycan separations. In addition, this column has the flexibility to be used in a purely ion-exchange, charge-based separation mode for the separation of various glycans without discrimination of size and polarity. This makes it a suitable tool for accurate quantification of glycans based on charge, which cannot be achieved with any other HPLC/UHPLC columns on the market. The GlycanPac AXH-1 column is designed for use with LC-fluorescence detection and LC/MS applications using volatile aqueous buffers (e.g., ammonium acetate or ammonium formate) and acetonitrile. The substrate of the GlycanPac AXH-1 column is based on high purity 1.9  $\mu\text{m}$  and 3.0  $\mu\text{m}$  spherical silica for UHPLC and standard HPLC applications, respectively.

In this application note we highlight the use of a GlycanPac AXH-1 column and Q Exactive hybrid quadrupole-Orbitrap mass spectrometer for the structural analysis of a 2-aminobenzamide (2AB) labeled *N*-glycans from bovine fetuin.

## Experimental Details

Consumables	Part Number
GlycanPac AXH-1, 1.9 $\mu\text{m}$ , analytical column (2.1 $\times$ 150 mm)	082472
Deionized (D.I.) water, 18.2 M $\Omega$ -cm resistivity	
Fisher Scientific™ HPLC grade acetonitrile (CH <sub>3</sub> CN)	AC610010040
Fisher Scientific LC/MS grade formic acid	A117-50
Fisher Scientific ammonium formate	AC40115-2500
Thermo Scientific Premium 2 mL vial convenience kit	60180-600
New England Biolabs® PNGase F	P0705L
Bovine fetuin from a commercial source	
Thermo Scientific™ Hypercarb™ 6 mL cartridge	60106-403
Fisher Scientific trifluoroacetic acid	28904
Fisher Scientific sodium cyanoborohydride	AC16855-0500
Fisher Scientific anthranilamide (2AB)	AC10490-5000
Fisher Scientific glacial acetic acid	AA36289AP

Equipment	Part Number
Thermo Scientific™ Dionex™ UltiMate™ 3000 Bio-RS system, consisting of LPG-3400RS pump, TCC-3000RS thermal compartment, WPS-3000TRS pulled-loop well plate autosampler, FLD3400RS fluorescence detector with Dual-PMT, and a 2 $\mu\text{L}$ micro flow cell	6078.4330
Q Exactive hybrid quadrupole-Orbitrap mass spectrometer	
Thermo Scientific™ SpeedVac™ system	
Thermo Scientific Lyophilizer (Labconco® FreeZone® -105 °C 4.5 L benchtop freeze dry system)	16-080-207
Thermo Scientific 24-port SPE vacuum manifold	60104-233

## Buffer Preparation

Ammonium formate (80 mM, pH 4.4): Dissolve 5.08  $\pm$  0.05 g of ammonium formate (crystal) and 0.60 g of formic acid in 999.6 g of D.I. water. Sonicate the resulting solution for 5 minutes.

## Sample Preparation

- Native *N*-glycans are released from glycoproteins with PNGase F enzyme and purified by a Hypercarb cartridge (6 mL) with the help of 24-port SPE vacuum manifolds under vacuum. The released glycans are conjugated with a 2-amino benzamide (2AB) label group using the reported procedure of Bigge et al. [7].
- Dissolve 2AB labeled *N*-glycan from fetuin (5,000 pmol) in 25  $\mu\text{L}$  D.I. water in a 250  $\mu\text{L}$  autosampler vial.
- Add 75  $\mu\text{L}$  acetonitrile to the same vial and mix till uniformity.

Note: Store the standard at -20 °C.

### Separation Conditions

Column:	GlycanPac AXH-1, 2.1 x 150 mm, 1.9 $\mu$ m					082472
Mobile phase A:	Acetonitrile / water (80:20, v/v)					
Mobile phase B:	Ammonium formate (80 mM, pH 4.4)					
Column temperature:	30 °C					
Sample volume:	1 $\mu$ L					
Gradient:	Time (min)	%A	%B	Flow Rate (mL/min)	Curve	
	-10	97.5	2.5	0.4	5	
	0	97.5	2.5	0.4	5	
	30	87.5	12.5	0.4	5	
	35	75.0	25.0	0.4	5	
	40	62.5	37.5	0.4	5	

### MS Conditions

MS instrument:	Q Exactive hybrid quadrupole-Orbitrap MS
Ionization mode:	Negative ion mode
MS scan range:	380-2000 $m/z$
Resolution:	70,000
AGC target	1 x 10 <sup>6</sup>
Max IT:	60 ms
dd-MS2 resolution:	17,500
MS/MS AGC target	2 x 10 <sup>5</sup>
MS/MS max IT:	1000 ms
Isolation window:	2 $m/z$
Dynamic exclusion:	90 s

### Data Processing and Software

Chromatographic software:	Thermo Scientific™ ChromQuest™ Chromatography Data System version 5.0
MS data acquisition:	Thermo Scientific™ Xcalibur™ software version 2.2 SP1.48
MS/MS data analysis:	SimGlycan® software (PREMIER Biosoft)

## Results

### Glycan Separation by Charge, Size and Polarity

Figure 1 shows the separation of neutral and acidic 2AB labeled *N*-glycans from bovine fetuin using a GlycanPac AXH-1 (1.9  $\mu$ m, 2.1 x 150 mm) column. The glycan elution profile consists of a series of peaks grouped into several clusters in which the neutral glycans elute first, followed by monosialylated, disialylated, trisialylated, tetrasialylated, and finally pentasialylated species. Analytes in each cluster represent glycans of the same charge. Within each cluster, the glycans having the same charge are further separated according to their sizes and polarity by HILIC interaction. The structure of the glycans present in each peak was determined in an LC-MS/MS study as shown in the following section.

### Structural elucidation

The 2AB labeled *N*-glycans from bovine fetuin were separated on the GlycanPac AXH-1 column based on the separation conditions using a two eluent system and analyzed on a Q Exactive benchtop mass spectrometer. The total ion chromatogram (TIC) is shown in Figure 1. For structural elucidation, data dependant MS/MS spectra were acquired on all precursor ions ( $z \leq 2$ ) and SimGlycan software from PREMIER Biosoft was used for data analysis [8]. The detailed structural information obtained (Table 1) from the MS/MS data further validated the ability of GlycanPac AXH-1 columns to separate glycans based on charge, size, isomers, and polarity. These results also confirmed that the GlycanPac AXH-1 column would be ideal for MS use.




Figure 1: LC/MS analysis of 2AB labeled *N*-glycans from bovine fetuin using the GlycanPac AXH-1 column

Peak	Compound Structure (2AB labeling are not shown)	Peak	Compound Structure (2AB labeling are not shown)
1		16	
2		17	
4		18	
5		19	
6		20	
7		21	
8		22	
9		23	
10		24	
11a		25	
11b		26	
11c			
12a			
12b			
13			
14			
15			

N-acetyl glucosamine (GlcNAc)

Mannose (Man)

Galactose (Gal)

N-acetyl neuraminic acid (Neu5Ac)

N-glycolyl neuraminic acid (Neu5Gc)

L-fucose (L-Fuc)

Table 1: Structural characterization of glycans present in each peak by the separation of 2AB labeled *N*-glycans from bovine fetuin using a GlycanPac AXH-1 column

## Conclusion

The GlycanPac AXH-1 column is a high-performance, silica-based HPLC column for simultaneous separation of glycans by charge, size, and polarity. It is designed for high-resolution and high-throughput analysis with unique selectivity for biologically important glycans. We have demonstrated that this column provides unique selectivity and excellent resolution for glycans released from fetuin.

LC-MS/MS analysis of 2AB labeled N-glycans derived from glycoproteins are demonstrated using GlycanPac AXH-1 columns. The Q Exactive hybrid quadrupole-Orbitrap instrument provides excellent MS/MS fragmentation information to enable characterization of glycan structures.

## Reference

- [1] Varki, A. Biological Roles of Oligosaccharides: All the theories are correct. *Glycobiology*, **1993**, *3*, 97–130.
- [2] Bertozzi, C. R., Freeze, H.H., Varki, A. and Esko, J.D. Glycans in Biotechnology and the Pharmaceutical Industry, Essentials of Glycobiology. 2nd edition. Chapter 51 Cold Spring Harbor (NY): Cold Spring Harbor Laboratory Press; 2009.
- [3] Guidance for Industry, Scientific Considerations in Demonstrating Biosimilarity to a Reference Product, 2012, U.S. ([www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM291128.pdf](http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM291128.pdf))
- [4] Thermo Scientific Accucore HPLC Columns Technical Guide: <http://thermo.dirxion.com/accucore>.
- [5] Thermo Scientific Accucore 150-Amide-HILIC application note on “Analysis of Human IgG Glycans on a Solid Core Amide HILIC Stationary Phase”: [http://www.separatedbyexperience.com/documents/AN-LC-Accucore-HILIC-Human-IgG-Glycans-AN20703\\_E.pdf](http://www.separatedbyexperience.com/documents/AN-LC-Accucore-HILIC-Human-IgG-Glycans-AN20703_E.pdf)
- [6] Thermo Scientific GlycanPac AXH-1 Product Specification: [http://www.dionex.com/en-us/webdocs/114170-PS-GlycanPac-AXH1-Column-PS20695\\_E.pdf](http://www.dionex.com/en-us/webdocs/114170-PS-GlycanPac-AXH1-Column-PS20695_E.pdf)
- [7] Bigge, J. C. et al., Non-selective and efficient fluorescent labeling of glycans using 2-amino benzamide and anthranilic acid. *Analytical Biochemistry*, **1995**, *230*, 229-238.
- [8] Apte, A and Meitei, N.S., Bioinformatics in glycomics: glycan characterization with mass spectrometric data using SimGlycan® software, *Meth. in Mol. Biol.*, **2010**, *600*, 269-81.



# Scantec Nordic

Analys & Mätteknik

031 336 90 00 • [www.scantecnordic.se](http://www.scantecnordic.se)

[thermoscientific.com/glycanpac](http://thermoscientific.com/glycanpac)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. New England Biolabs is a registered trademark of New England Biolabs. Labconco and FreeZone are registered trademarks of Labconco Corp. SimGlycan is a registered trademark of PREMIER Biosoft. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**United Kingdom** +44 (0)1928 534110  
**Japan** +81 3 5826 1615

**China** +86 21 68654588 +86 10 84193588  
 +86 20 83145199 800 810 5118  
**India** +91 22 6742 9494 +91 27 1766 2352  
**Australia** 1 300 735 292 (free call domestic)  
**New Zealand** 0800 933 966 (free call domestic)  
**All Other Enquiries** +44 (0) 1928 534 050

**Technical Support**  
**North America** +1 800 332 3331  
**Outside North America** +44 (0) 1928 534 440

**Thermo**  
 S C I E N T I F I C

Part of Thermo Fisher Scientific