

Development of Ultra-fast pH-Gradient Ion Exchange Chromatography for the Separation of Monoclonal Antibody Charge Variants

Ken Cook,¹ Frank Steiner,² Mauro De Pra²

¹Thermo Fisher Scientific, Hemel Hempstead, United Kingdom; ²Thermo Fisher Scientific, Germering, Germany

Key Words

Critical Quality Attributes, Biotherapeutics, Intact Proteins, Vanquish UHPLC, MAbPac SCX-10 RS Column

Introduction

The pre-formulated buffers for pH gradient, introduced by Thermo Fisher Scientific, have greatly simplified the development of ion exchange chromatography (IEX) of monoclonal antibodies (mAbs). Three features make this simplification possible. The first feature is that the buffers can cover a pH range from 5.6 to 10.2; this pH window enables the characterization of mAbs with a wide range of isoelectric points, providing a global pH gradient ion exchange screening method that will accommodate the majority of therapeutic mAbs. The second feature is that, if a linear solvent gradient is programmed in the pump, the actual pH gradient produced in the column will be linear as well. Although this sounds trivial, scientists involved in pH gradient studies know how difficult it is to develop buffer formulations capable of fulfilling this requirement. The advantage of a genuine linear pH gradient is that the method can be confidently fine-tuned merely by narrowing down the pH range around the mAb and its variants, thus allowing for the adjustment of the gradient time according to the resolution requirements. The third feature is the mobile phase preparation: the pre-formulated pH buffers only need to be diluted by a factor of 10 in deionized water, and the mobile phase is ready to use. Compared to a salt gradient ion exchange chromatography method, where no generic screening can be easily designed, and where method development goes through the rather tedious preparation of several buffers at different pH values, the time and effort invested in method development are substantially reduced.



High resolution pH gradient separations are obtained with 30 minute gradients and relatively long columns, such as the Thermo Scientific™ MAbPac™ SCX-10 column, 10 μ m, 4 \times 250 mm.¹ Fast separation of mAb variants were demonstrated by using the MAbPac SCX-10 column, 5 μ m, 4 \times 50 mm with no significant loss in resolution² on a Thermo Scientific™ UltiMate™ 3000 BioRS system. In this case, the pH gradient from 5.6 to 10.2 was completed in 7.5 minutes, and the resolution between variants was still satisfactory, despite the short analysis time.

In this work, we show how to push the throughput of pH gradient IEX even further. To achieve this, a MAbPac SCX-10 RS column, 5 μ m, 2.1 \times 50 mm was operated on a new Thermo Scientific™ Vanquish™ UHPLC system. The Vanquish UHPLC is a fully biocompatible system, suitable for the analysis of intact proteins. The combination of low gradient delay volume and high precision gradient formation makes it the ideal system for high throughput analysis with gradient elution. Here the system was used with the default configuration and a total system gradient delay volume of 175 μ L. Fast charge variant separations of 5 mAbs are shown.

Goal

Provide high throughput pH gradient separations of mAb charge variants

Experimental

Vanquish UHPLC, consisting of:

- Vanquish System Base (P/N VH-S01-A)
- Binary Pump H with Default Mixer (P/N VH-P10-A)
- Split Sampler HT (P/N VH-A10-A)
- Column Compartment H (P/N VH-C10-A)
- Diode Array Detector HL (P/N VH-D10-A)

Chromatographic Conditions

| | |
|--------------------------------|---|
| Column: | MABPac SCX-10 RS, 5 μ m, 2.1 \times 50 mm (P/N 082675) |
| Buffers: | Thermo Scientific CX-1 pH-Gradient buffer A (pH 5.6) 125 mL (P/N 083273) CX-1 pH-Gradient buffer B (pH 10.2) 125 mL (P/N 083275) |
| Mobile Phase A: | CX-1 pH-Gradient buffer A (pH 5.6) diluted 10x in deionized water |
| Mobile Phase B: | CX-1 pH-Gradient buffer B (pH 10.2) diluted 10x in deionized water |
| Column Compartment Temperature | 30 $^{\circ}$ C, forced air |

Detector and Conditions

| | |
|-------------------------|---|
| Detector: | LightPipe™ 10 mm Standard Flowcell (P/N 6083.0100) |
| Detection Wavelength: | 280 nm |
| Data Acquisition Range: | 5 Hz (for flow rate \leq 0.5 mL/min) and 50 Hz (for flow rate \geq 1.0 mL/min) |
| Response Time: | 2 s (for flow rates \leq 0.5 mL/min) and 0.1 s (for flow rates \geq 1.0 mL/min) |

Data Processing

| | |
|-----------|---|
| Software: | Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Chromatography Data System (CDS) |
|-----------|---|

Samples

| mAb Name | Concentration (mg/mL) | Injection Volume (μ L) |
|-------------|-----------------------|-----------------------------|
| Bevacizumab | 25 | 1 |
| Cetuximab | 5 | 4 |
| Infliximab | 10 | 4 |
| Trastuzumab | 21 | 2 |
| mAb A | 21 | 2 |



Results and Discussion

One of the benefits of using CX1 buffers for the pH-gradient is the simplified method development and optimization. This is due to the fact that the pump running a linear solvent gradient will result in a linear pH gradient in the column. This is not the case for most of home-made buffer formulations, which would produce a non-linear pH gradient in response to a linear programmed gradient. A non-linear pH gradient makes method optimization difficult due to the uncertainty of the actual effects of any changes in the programmed gradient.

It is recommended to perform a generic screening from pH 5.6 to 10.2 when the pH at which a given mAb elutes is not known. In this work, the generic screening was run in 10 minutes at 0.45 mL/min. As it can be seen in Figure 1 in some cases, satisfactory separation of the charge variants was achieved during the first run. This was the case for cetuximab and infliximab: several charge variants could be resolved with sufficient resolution.




Figure 1. Separation of 5 mAbs with a generic 0 \rightarrow 100 %B in 10 min at 0.45 mL/min (method #1 Table 1). Samples are as follows: A) bevacizumab, B) cetuximab, C) infliximab, D) trastuzumab, E) mAb A.

Table 1. Description conditions used for the separation of charge variants of the different monoclonal antibodies.

| Method # | Figure | Gradient Range (%B) | Gradient pH Range | Gradient Time (min) | Flow Rate (mL/min) | Gradient Slope* (%B/mL) |
|----------|--------|---------------------|-------------------|---------------------|--------------------|-------------------------|
| 1 | 1 | 0–100 | 5.6–10.2 | 10.0 | 0.45 | 22.2 |
| 2 | 2A | 20–40 | 6.5–7.4 | 5.0 | 0.50 | 8.0 |
| 3 | 2B | 23–35 | 6.7–7.2 | 2.5 | 1.00 | 4.8 |
| 4 | 3A | 10–35 | 6.1–7.2 | 5.0 | 0.50 | 10.0 |
| 5 | 3B | 10–35 | 6.1–7.2 | 2.5 | 1.00 | 10.0 |
| 6 | 4A | 20–40 | 6.5–7.4 | 5.0 | 0.50 | 8.0 |
| 7 | 4B | 18–27 | 6.4–6.8 | 0.8 | 1.20 | 9.4 |
| 8 | 5A | 35–60 | 7.2–8.4 | 5.0 | 0.50 | 10.0 |
| 9 | 5B | 33–45 | 7.2–7.7 | 2.5 | 1.00 | 4.8 |
| 10 | 6B | 5–30 | 5.8–7.0 | 2.5 | 1.00 | 10.0 |

* Slope based on gradient volume

After the generic screening, the method development efforts aimed at decreasing the analysis cycle time and at the same time improving resolution. To achieve this goal, two parameters were modified, namely the pH range and the gradient slope. Here we used gradient slope based on gradient volume, i.e. $\Delta(\%B)/V_G$, where %B is the amount of B eluent and V_G is the volume of mobile phase delivered by the pump during the gradient. A narrower pH window allowed for a reduced run-time, whereas a shallower gradient slope provided better resolution.

The gradient slope of the generic screening between pH 5.6 and 10.2 was 22.2 (%B)/mL; the improved and faster analysis were obtained with gradient slopes of 8 (%B)/mL and 10 (%B)/mL. Details of the conditions are listed in Table 1. Figures 2a, 3a, 4a, and 5a, show the results obtained by this approach. The number of resolved variants was larger than in the initial screening. Gradient time was 5 minutes, hence half of the time of the initial screening.




Figure 2. Separation of bevacizumab with different gradient conditions. Chromatogram A: method #2 Table 1; chromatogram B: method #3 Table 1.




Figure 3. Separation of cetuximab with different gradient conditions. Chromatogram A: method #4 Table 1; chromatogram B: method #5 Table 1.




Figure 4. Separation of infliximab with different gradient conditions. Chromatogram A: method #6 Table 1; chromatogram B: method #7 Table 1.

The following step was used to develop fast analysis cycles compatible with high throughput. With this approach we aimed for a 2.5 minute gradient time or lower, and total analysis time of less than 4 minutes, including column re-equilibration. Data are shown in Figures 2b, 3b, 4b, 5b, and 6b. The purpose was to develop a method suitable for high-throughput that can run at least 300 samples a day. High flow rate was used for this purpose. A high flow rate allows running short gradients with relatively shallow gradient slopes; the shallow slope is required to preserve selectivity between charge variants. Additionally, column equilibration, which is directly dependent on the volume of mobile phase flowing through the column, is reached quicker. MAbPac SCX-10 RS column is pressure rated up to 7000 psi (~ 480 bar), therefore it can be operated at high linear flow rate. In this work, we used flow rates up to 1.2 mL/min.

The chromatographic pattern between methods at moderate and high flow rate was preserved. The separation capabilities of different methods were compared based on the resolution between charge variants. Since in several instances peak pairs were overlapping, and it was not always possible to measure peak width at half height or at the baseline, here we used resolution based on statistical moments. Resolution was calculated directly by Chromeleon 7.2 CDS according to the formula:

$$R = \frac{(t_R^2 - t_R^1)}{2 * (\sqrt{\mu_2^2} + \sqrt{\mu_1^2})}$$

where t_R^2 and t_R^1 are the retention times of the more and less retained peak respectively, and μ_2^2 and μ_1^2 are the related second moment.

In some cases, the ultra-fast separation approach was accompanied by some resolution loss. This is the case of trastuzumab, where average resolution loss was ~ 3%. In the case of the complex variants pattern of cetuximab, the average resolution loss at high flow rate was ~ 13%. Infliximab resolution decreased by ~ 11%, however separation of the 5 main charge variants and 2 minor ones was achieved in 1 minute. This impressive result was obtained by running the column at 1.2 mL/min with a 0.8 min gradient time. In the case of bevacizumab, the ultra-fast separation approach even yielded to ~ 4% improved resolution. The better resolving power can be explained by a slightly narrower pH range and a shallower gradient slope.

mAb A was analyzed only with the generic screening and high flow rate (1 mL/min) method. The high-throughput method provided 19% better resolution on average of 5 charge variants.




Figure 5. Separation of trastuzumab with different gradient conditions. Chromatogram A: method #8 Table 1; chromatogram B: method#9 Table 1.




Figure 6. Separation of mAb A with different gradient conditions. Chromatogram A: method #1 Table 1; chromatogram B: method#10 Table 1.

Table 2. Overview of the resolution between charge variants at different conditions. Resolution was calculated by Chromeleon CDS using statistical moments.

| Statistical Resolution Charge Variants | | | | | | | |
|--|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|
| mAb | Gradient Time (min) | Peaks 1–2 | Peaks 2–3 | Peaks 3–4 | Peaks 4–5 | Peaks 5–6 | Peaks 6–7 |
| Bevacizumab | 5 | 0.94 | 0.87 | 1.81 | — | — | — |
| | 2.5 | 1.00 | 0.89 | 1.89 | — | — | — |
| Cetuximab | 5 | 0.79 | 1.14 | 1.00 | 1.13 | 1.28 | 0.98 |
| | 2.5 | 0.53 | 1.05 | 0.94 | 0.93 | 1.20 | 0.91 |
| Infliximab | 5 | 0.68 | 1.56 | 1.33 | 0.87 | — | — |
| | 0.8 | 0.59 | 1.40 | 1.30 | 0.70 | — | — |
| Trastuzumab | 5 | 1.25 | 1.14 | 1.70 | 0.87 | — | — |
| | 2.5 | 1.27 | 1.11 | 1.63 | 0.79 | — | — |
| mAb A | 10 | 0.14 | 1.31 | 0.94 | 1.26 | — | — |
| | 2.5 | 0.32 | 0.96 | 1.22 | 2.02 | — | — |

Conclusion

The ultra-fast charged variant separations described here are achieved because of several advances in chromatography techniques. The mechanism of pH gradient chromatography lends itself to the use of shorter, faster columns. The availability of high pressure rated small particle size ion exchange columns are a perfect match to pH gradient methodology. The commercial buffer formulations used here form a linear gradient which allows intelligent optimization of the methods. Finally, there is the use of the new Vanquish UHPLC system which has extremely low delay volumes, high precision gradient formation and a totally inert flow path.

Ultra-fast separation that requires total analysis cycle in the order of 2 minutes, including column re-equilibration and injection time, enables users to run more than 1,400 samples during 48 hours continuous operations. To allow unattended tasks with such large amount of samples, the Vanquish UHPLC system can be extended with the Vanquish Charger Module. This can host up to 9000 samples in a thermostatted environment, and transfer them to the Vanquish Autosampler.

References

1. Thermo Fisher Scientific. Thermo Scientific Application Note 20784, *A Novel pH Gradient Separation Platform for Monoclonal Antibody (MAb) Charge Variant Analysis*, AN20784_E, Sunnyvale, CA, 2013.
2. Thermo Fisher Scientific. Thermo Scientific Application Note 20946, *A Fast and Robust Linear pH Gradient Separation Platform for Monoclonal Antibody (mAb) Charge Variant Analysis*, AN20946_E, Sunnyvale, CA, 2014.

www.thermoscientific.com/chromatography

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific,
Sunnyvale, CA USA is
ISO 9001:2008 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand